
Sequencing by synthesis



FASTQ files

Line1: Sequence identifier


Line2: Raw sequence

Line3: meaningless

Line4: quality values for the sequence

```
@HWI-ST508:210:C0EDTACXX:1:1101:1872:1227 1:N:0:
AATTGTGAAAACCCAAAAGGTGGAGCAGCCATTNTTATACATTGCAGAAGGGNGANNNANCNTTATGAAATTTAGCACCTGCCTTCCTGAATGATAAATGG
@CCFFEFFHHHHHJJJJIJJCGHEIIIJIJJJJ#1BFHIJJJJJJJJJJJJJI#-;###-#-#-5?BFFFFEEEEEECCDDDDDDDDDDDCCDDDDDCCEED
@HWI-ST508:210:COEDTACXX:1:1101:1895:1233 1:N:0:
IGACATAAGCTTGCATTTGAAAAGCACCTCCGAAAGCTTCCCAGCCTCAAAGNCANNATCGNCTTCTGATGCAGTTAGGCACCACAAGAGCTTCCCCACAA
@HWI-ST508:210:C0EDTACXX:1:1101:1761:1235 1:N:0:
GCTCTACTAAAAATATAAAAATTGGCCAGGCGCAGTGACACATGCCTGTAGTCCCNGCTATTCGGGAGGCTGACACACAAGAATCAATCACTTGAACCCAG
CCCFFFFFHGHHHJJJJJJJJJJJJJJJJIEIIIJFHGIIIIJJIJJJJHIJJIJ#-;FGGIJIJHHFFDDEEDDCCDDDDCCDDDDDDDDDDDDDDDD
@HWI-ST508:210:C0EDTACXX:1:1101:1971:1236 1:N:0:
CAGGATGAAAGAGGTCTGGCCAGGTGCTGGGTGCAGTGGCTCACACCTGTAATCCCAGCACTTTGGGAGGCCGAGGTGGGCGGATCACGAAGTCAGGAGTT
CCCFFFFFHGHHGJHIJIIJJJJI3CFGIJJ9DFHJDEHGIJIJJJJJIIJJJGGIJJJJJJJFIJHFFFFDDDB/?BB@BD<39?CD@B8+:@CDCB##
@HWI-ST508:210:C0EDTACXX:1:1101:1830:1239 1:N:0:
@HWI-ST508:210:C0EDTACXX:1:1101:1999:1240 1:N:0:
@0@0 	ext{DDA2} FHBHHEGEHIHGIGGHBFCGIEHGAEGGIIEGIIIIGHIGEHEGHIGIGBFHEHIEAHGHHFHEH; B@0 	ext{EBDCDEEBCDDCCCCC}
@HWI-ST508:210:C0EDTACXX:1:1101:1806:1245 1:N:0:
```

FastQC



http://www.bioinformatics.babraham.ac.uk/projects/fastqc/

SAM format


• Standard format for short-read alignements.

```
HWI-ST151_106137860:1:67:20248:73945\#0 129 chr17 98508 255 40M =
```


98849 378 AGGGGTTGGCGGGGCAAGGTGGCTCACGCCTGTCATCCCA

@B@B@:8A?8>@@80DCCDA@85,C7>7>>AB#########

Visualization



Visualization



VCF format

```
#fileformat=VCFv4.0
##FILTER=<ID=ABFilter,Description="AB > 0.75">
##FILTER=<ID=HRunFilter,Description="HRun > 5.0">
##FILTER=<ID=LowQual, Description="Low quality">
##FILTER=<ID=ODFilter, Description="OD < 5.0">
##FILTER=<ID=OUALFilter,Description="OUAL < 30.0">
##FILTER=<ID=SBFilter,Description="SB > -0.10">
##FILTER=<ID=SnpCluster.Description="SNPs found in clusters">
##FORMAT = < ID = AD, Number = ., Type = Integer, Description = "Allelic depths for the ref and alt alleles in the order listed" >
##FORMAT = <ID=DP, Number=1, Type=Integer, Description="Read Depth (only filtered reads used for calling)">
##FORMAT=<ID=GO.Number=1.Type=Float.Description="Genotype Quality">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=PL, Number=3, Type=Float, Description="Normalized, Phred-scaled likelihoods for AA, AB, BB genotypes where A=ref and B=alt; not applicable
##INFO=<ID=AB, Number=1, Type=Float, Description="Allele Balance for hets (ref/(ref+alt))">
##INFO=<ID=AC, Number=., Type=Integer, Description="Allele count in genotypes, for each ALT allele, in the same order as listed">
##INFO=<ID=AF, Number=., Type=Float, Description="Allele Frequency, for each ALT allele, in the same order as listed">
##INFO=<ID=AN, Number=1, Type=Integer, Description="Total number of alleles in called genotypes">
##INFO=<ID=BaseQRankSum,Number=1,Type=Float,Description="Phred-scaled p-value From Wilcoxon Rank Sum Test of Alt Vs. Ref base qualities">
##INFO=<ID=DB, Number=0, Type=Flag, Description="dbSNP Membership">
##INFO=<ID=DS, Number=0, Type=Flag, Description="Were any of the samples downsampled?">
##INFO=<ID=Dels, Number=1, Type=Float, Description="Fraction of Reads Containing Spanning Deletions">
##INFO=<ID=HRun, Number=1, Type=Integer, Description="Largest Contiguous Homopolymer Run of Variant Allele In Either Direction">
##INFO=<ID=HaplotypeScore.Number=1.Type=Float.Description="Consistency of the site with at most two segregating haplotypes">
##INFO=<ID=MQ, Number=1, Type=Float, Description="RMS Mapping Quality">
##INFO=<ID=MO0, Number=1, Type=Integer, Description="Total Mapping Quality Zero Reads">
##INFO=<ID=MQRankSum, Number=1, Type=Float, Description="Phred-scaled p-value From Wilcoxon Rank Sum Test of Alt Vs. Ref read mapping qualities">
##INFO=<ID=OD, Number=1, Type=Float, Description="Variant Confidence/Quality by Depth">
##INFO=<ID=ReadPosRankSum, Number=1, Type=Float, Description="Phred-scaled p-value From Wilcoxon Rank Sum Test of Alt Vs. Ref read position bias">
##INFO=<ID=SB, Number=1, Type=Float, Description="Strand Bias">
##UnifiedGenotyper="analysis type=UnifiedGenotyper input file=[bamFiles.list] sample metadata=[] read buffer size=null phone home=STANDARD read f
##VariantFiltration="analysis type=VariantFiltration input file=[] sample metadata=[] read buffer size=null phone home=STANDARD read filter=[] in
##GATK version=1.0.5614
##Reference Sequence=<url=https://grc-aspera-public:Asppass1@aspera.gs.washington.edu/aspera/user/?B=%2Fhuman refseq, file=hg19 genome reference.
##DBSNP ROD=<url=https://grc-aspera-public:Asppass1@aspera.gs.washington.edu/aspera/user/?B=%2Fdbsnp, file=chr order.sorted.dbsnp 131 hg19.rod.gz
##Exome Target=<url=https://grc-aspera-public:Asppass1@aspera.gs.washington.edu/aspera/user/?B=%2Fexome targets, file=nimblegen solution uwrefseg
#CHROM POS
               ID
                                        OUAL FILTER INFO FORMAT 1796 1797 1798 1799 1800 1801 1802 1803 1804
                                        752.59 QDFilter;SBFilter
        69270 .
                                                                        AC=32:AF=0.667;AN=48;BaseORankSum=53.792;Dels=0.00;HRun=0;HaplotypeScore=
        69428 rs71245814
                                                2215.91 PASS AB=0.61:AC=4:AF=0.0263:AN=152:BaseORankSum=200.000:DB:DS:Dels=0.00:HRun=0:Haplotv
        69511 rs2691305
                                                43693.52
                                                                PASS AB=0.61;AC=110;AF=0.7857;AN=140;BaseORankSum=19.099;DB;DS;Dels=0.00;HRun=
                                        406.53 SBFilter
                                                                AC=2; AF=0.0149; AN=134; BaseQRankSum=50.877; Dels=0.00; HRun=0; HaplotypeScore=0.3804;
        69680 .
        69897 rs75758884
                                                340.36 SBFilter
                                                                         AC=22; AF=0.846; AN=26; BaseQRankSum=21.543; DB; Dels=0.00; HRun=1; HaplotypeSco
                                        245.91 PASS
        865584 .
                                                        AB=0.56; AC=2; AF=0.0105; AN=190; BaseORankSum=9.661; Dels=0.00; HRun=0; HaplotypeScore=0.3829; M
        865628 rs41285790
                                                383.75 PASS AB=0.57; AC=1; AF=0.0053; AN=190; BaseQRankSum=43.179; DB; Dels=0.00; HRun=0; HaplotypeSc
                                                                AB=0.53; AC=14; AF=0.0737; AN=190; BaseQRankSum=200.000; DB; Dels=0.00; HRun=0; Haplotype
        865694 rs9988179
                                                7122.84 PASS
        865700 . C
                                        646.66 PASS
                                                        AB=0.50; AC=2; AF=0.0105; AN=190; BaseORankSum=9.945; Dels=0.00; HRun=0; HaplotypeScore=1.3140; M
        866422 .
                                        1045.98 PASS
                                                        AB=0.59; AC=1; AF=0.0053; AN=190; BaseQRankSum=200.000; DS; Dels=0.00; HRun=1; HaplotypeScore=3.7
        866438 .
                                                        AB=0.48; AC=1; AF=0.0053; AN=190; BaseQRankSum=3.012; DS; Dels=0.00; HRun=0; HaplotypeScore=2.887
                                        2011.24 PASS
```

VCF

```
##fileformat=VCFv4.2
##contig=<ID=2,length=51304566>
##INFO=<ID=AC, Number=A, Type=Integer, Description="Allele count in genotypes">
##INFO=<ID=AN, Number=1, Type=Integer, Description="Total number of alleles in called genotypes">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=DP, Number=1, Type=Integer, Description="Read Depth">
##FORMAT=<ID=GQ, Number=1, Type=Integer, Description="Genotype Quality">
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT
                                                                             SAMPLE3
                                                                                         SAMPLE4
                                                                                                     SAMPLE5
                                                                                                                 SAMPLE6
                                                                                                                             SAMPLE7
                                                     SAMPLE1
                                                                 SAMPLE2
2 81170 . C T
                                          GT:DP:GQ
                                                    0/0:4:12
                                                                 0/0:3:9
                                                                             0/1:1:3
                                                                                         0/1:9:24
                                                                                                     1/0:4:12
                                                                                                                 0/0:5:15
                                                                                                                             0/0:4:12
                          AC=9; AN=7424
                                                    0/1:4:12
                                                                             0/0:1:3
  81171 . G A
                          AC=6; AN=7446
                                          GT:DP:GQ
                                                                 0/0:3:9
                                                                                         0/0:9:24
                                                                                                     0/1:4:12
                                                                                                                 0/1:5:15
                                                                                                                             0/0:4:12
                          AC=5; AN=7506
  81182
          . A G
                                          GT:DP:GQ
                                                    0/0:5:15
                                                                 0/0:4:12
                                                                             0/0:5:15
                                                                                         0/0:9:24
                                                                                                     0/0:4:12
                                                                                                                 0/0:4:12
                                                                                                                             0/0:4:12
  81204 . T G
                          AC=2; AN=7542
                                          GT:DP:GQ
                                                    1/0:5:15
                                                                 0/0:9:27
                                                                             0/0:10:30
                                                                                         0/0:15:39
                                                                                                     0/0:9:27
                                                                                                                 1/0:13:39
                                                                                                                             0/1:14:42
                  . .
```

Primary, secondary and tertiary analysis

Primary analysis

Happens in the instrument, checks during the run for quality control. Consists of image alignment, color and quality value calls. **FASTQ files**

Secondary analysis

A reference sequence is converted into color space and the data are aligned to the reference sequence.

A consensus sequence may then be constructed from the sequencing reads.

Comparison of the consensus sequence to a reference genome enables the identification of SNPs and structural variations. SAM/BAM files

<u>Tertiary analysis</u>

Data analysis that takes place after the reads are mapped. Visualization of data. Annotation of variants. VCF files